



Whole You Sleep Appliance Respire Pink Micro




Shorter Arms for Increased Wearable Comfort Even for Patients with a Small Mouth

Three reasons to choose the Respire Pink Micro

1 Shorter Adjustment Arms



2 Freedom to Move Your Mouth



3 Easily Adjustable



Respire Pink Micro Details

- Durability suitable for patients who grind their teeth at night
- 6 mm maximum advancement with an extra anterior fixing element
- Freedom in mouth movement allows patients to move their lower jaw vertically and laterally
- Available in hard and hard/soft-fitting surface
- 1-year manufacturer's warranty
- Option to add elastic hooks and anterior disclosing ramp
- Easily adjustable
- FDA cleared